

Apollo

Sitting Stool

 statmedical
we take care



Enhance Your Workspace with the Apollo Pneumatic Stool: Superior Medical-Grade Seating

Discover the Apollo Pneumatic Stool, our leading-edge release designed with Canadian hospitals and their rigorous infection control standards in mind. With the Apollo, you are not getting a simple office chair but a powerful commitment to safety, comfort, and efficiency.



Features

- Medical-grade seating upholstered with industry-best materials that offer unparalleled durability and resistance against chemicals
- High-risk areas are safeguarded with robust plastic to avert potential damage
- Adjustable back and height and removable arms for maximum ergonomic comfort

Options

Customizable details include:

- Toe Cap base style
- Stationary glide
- Soft Wheel Casters
- Auto-Lock Casters
- Reverse Auto-Lock casters

	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
Apollo-LAB	27"	27"	18.5"	41-51"	23-33"	300 lbs



Apollo



Apollo-AL



Apollo-P