

# Cama Recliner

 statmedical  
we take care



# Independent Reclining for Therapeutic Comfort

The smooth reclining motions, mixed with the iconic Cama stylized therapeutic comfort, make the Cama recliner an excellent addition to the patient's recovery room. Foot and backrest controls are independent, allowing for infinite positions, from fully upright to a flat sleeping setting. Its composite foam seat with dynamic suspension ensures long-term comfort, and a variety of options make it a versatile addition.



## Features

- Independent back and footrest allow for infinite reclining positions, including lay-flat
- Composite foam seat with a dynamic suspension system for long-term comfort
- Polyester caps resistant to hospital cleaners
- Crumb space in the seat for an easier cleanout

## Options

- Drainage bag hook-up that comes in a set of two
- Oxygen tank holder and IV Pole
- Seat moisture barrier to enhance durability
- Cal 133 flammability
- Central Locking Casters
- Folding arm and rolling foot tray
- Ergotron TV arm/no wiring
- Heated seat & lumbar or lumbar with massage
- Continuous footrest for added comfort
- BKH Swivel profile tablet - thermofoil top
- Folding/Swiveling side tablet - thermofoil top
- Power outlets and USB port options

	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
108-56	37"	32"	24"	45"	18"	550 lbs

