

# Eclipse

## Overbed Table

 statmedical  
we take care



# Supporting Wellness with Smart Design and Accessible Cost

Elevate your space with a practical, budget-friendly overbed table that doesn't sacrifice performance. Introducing the Eclipse Overbed Table, a cost-effective solution with high quality. Engineered to meet tight budgets, this table boasts an attractive thermoformed top designed to contain large spills. The sturdy powder-coated H-base ensures stability, while the spring-assisted lift mechanism provides infinite height adjustment from 29" to 45".



## Features

- Spring-assisted lift mechanism with infinite height positioning
- Thermoformed top with spill containment surface
- Powder-coated steel column and base for durability

## Options:

- Vanity tray with mirror
- Wire storage basket
- Choose from various powder-coat or chrome finishes to match your aesthetic

	Depth	Width	Height	Weight Capacity
ST140	16"	30"	29-45"	25 lbs

