

Flo BH

Room Furniture

 statmedical
we take care



Unmatched Durability and Functionality for Behavioral Healthcare Spaces

Flo BH Series is a rotationally molded case goods collection meticulously designed for behavioral healthcare environments. Comprising a bed, bedside units, and shelving, each piece is crafted with high-impact polyethylene, ensuring unparalleled durability through continuous molding. Flo not only withstands the rigors of behavioral health settings but also provides a seamless blend of functionality and longevity.



Flo BH

Room Furniture



Features

- The bedside unit incorporates integrated rails for secure containment on the top surface
- Bed and cabinet sides fit flush to walls, optimizing space utilization
- Bed features surface grooves for efficient fluid management
- Single-side bed storage can be positioned left or right-facing
- Stain-resistant, seamless, and easy-to-clean surfaces
- Corners and edges are rounded for increased safety
- 100% recyclable material



Options

- Enclosed bottom for added versatility
- Secure units in place with bolt-to-floor brackets for stability
- Patient restraint module for safety
- Non-removable nylon glides for chairs and tables
- More Bed options available to fit your requirements:
 - Bolt-to Floor
 - Enclosed Bottom
 - Patient Restraint Module
 - Corner Wall Filler
 - Head Wall Filler
 - Mattress
 - Mattress with Pillow

FLBD3680

Depth 83.5"
Width 39.5"
Height 20"
Top Height 18"
Weight 146 lbs



FRCG702

Depth 14"
Width 29"
Height 22"
Top Height -
Weight 23 lbs



FLBDS3680

Depth 83.5"
Width 39.5"
Height 20"
Top Height 18"
Weight 155 lbs



FRCG703

Depth 14"
Width 29"
Height 51"
Top Height -
Weight 49 lbs



FL100

Depth 18"
Width 20"
Height 23.4"
Top Height 22.4"
Weight 26 lbs



SF36805

Depth 36"
Width 80"
Height 5"
Top Height -
Weight -



SF36805-PW

Depth 36"
Width 80"
Height 5"
Top Height -
Weight -



