

Oasis

Treatment Recliner

 statmedical
we take care



Elevate Care with Oasis Treatment Recliner

The Oasis Treatment Recliner is a haven of comfort meticulously designed for the demands of 24/7 healthcare environments. This recliner, true to its name, offers a serene retreat with generous padding and an infinite recline feature. In the often non-relaxing atmosphere of hospitals, the Oasis Treatment Recliner stands out as a comforting refuge.



Features

- Easily maintain and refresh with field-replaceable upholstery covers
- Empower caregivers and patients with a side lever for footrest operation
- Infinite-position recline mechanism
- Smooth and secure transport with ergonomic handles

Options

- Available with upholstered arms, polyurethane, or wood arm caps
- Therapeutic benefits with Trendelenburg positioning
- Added convenience with a slide-out foot tray
- Opt for a moisture barrier for increased durability
- Extra comfort with upholstered arms
- Customize functionality with an IV pole and side tables

	Depth	Width	Seat Width	Height	Seat Height	Weight capacity
SR900-23	36"	29"	23"	48"	20.5"	500 lbs
SR900-26	36"	32"	26"	48"	20.5"	500 lbs



SR900-23



SR900-26