

Olsen

Overbed Table

 statmedical
we take care



Enhance Patient Experience with Olsen Overbed Table

Introducing the Olsen Overbed Table, a visually appealing and incredibly functional addition to healthcare spaces. With seamless and quiet height adjustability ranging from 28" to 45", Olsen ensures a customized and comfortable experience for patients. The optional flip-up extension, providing an extra 10 inches of surface area, adds to its versatility. Choose Olsen for a perfect blend of style and practicality.



Features

- Pneumatic cylinder mechanism for infinite height positioning
- Low profile U-base and C-base compatible with various medical equipment
- Dual-wheel casters for quiet maneuverability
- Spill containment top surface for added convenience

Options

- Convenient flip-up extension for added surface space
- Wire basket for easy storage and accessibility
- Vanity mirror for patient convenience
- Choice of thermofoil or solid surface top

	Depth	Width	Height	Base Height
ST200	32"	18"	28" – 45"	2.125"
ST200-EXT	44"	18"	28" – 45"	2.125"



ST200



ST200-SO



ST200-EXT



ST200-EXT-SO