

Split-Top Overbed Table

 statmedical
we take care



Extended Top, Increased Function

The Split-Top Overbed Table provides extra functionality and is designed with a secondary top that slides to either side of the primary top. A tab on the structure releases the secondary top, which is so light it can be moved with a single hand. This table maintains balance, regardless of where the secondary top is moved to, securing a better experience for patients and caregivers alike.



Features

- Thermoformed top with spill containment
- Pneumatic cylinder
- Infinite height adjustment between 29.5" and 44.5"
- Oversized 2" dual-wheel casters
- Slide-out top for extra surface area

Options

- Single-Wheel Rubber Casters
- H- or U- style base
- Column on left or right side
- Vanity Tray

Depth

Width

Height

Weight Capacity

ST170

18"

32"

29.5" - 44.5"

25 lbs

