

Transcend Overbed Table

 **statmedical**
we take care



Sleek, Durable, and Versatile Table for Healthcare Spaces

Find the perfect blend of functionality and modern design with the Transcend Overbed Table. Its sleek lines, rugged durability, and innovative features make it an ideal choice for healthcare environments. The low-profile U-base easily fits under late-model beds, stretchers, mobile recliners, and wheelchairs, while the pneumatic cylinder ensures smooth and nearly silent height adjustments.



Features

- Gas cylinder height adjustment mechanism for precise control
- Impressive 25 lb weight capacity for versatile use
- Variable height positioning ranging from 29.5" to 44.5"
- Anodized aluminum column for durability
- 14-gauge steel base covered with molded ABS shroud for protection

Options

- Thermoform or poly-resin top available
- Vanity tray with a mirror for added convenience
- Custom poly-resin edge colors with a minimum 30-unit order

	Depth	Width	Height	Weight Capacity
ST180	18"	32"	29.5" - 45.5"	25 lbs

